

Procedure Steps


Opening Software / Prepare Hardware

- Oculus App (Shortcut on Desktop)
 - Only needs to be opened - No need to press or do anything additional here
 - Having this application open allows us to enable link connections later
- Unity Hub (Shortcut on Desktop)
 - Open the most recent Unity project titled “assemblyTask - unity”
 - There are two with the same name located here
 - Confirm the file location as needed:
“C:\Users\shine\Development\git\Jalynnn\asTask”
 - This will then open the corresponding Unity Application
- Unity Application (Opened through Unity Hub)
 - Confirm that you are in the correct scene: “Assets > Scenes > New Scenes > Tutorial Video”
 - You can also confirm by checking the far top left of the monitor to confirm that it says “Tutorial Video”
 - Do **not** yet press play
- Locate LiveAmp EEG Device
 - Turn on the device by holding the power button until the blinking green turns solid green.
 - Enter pairing mode; a light under the green will turn solid blue
- Open LiveAmp.exe (Shortcut on Desktop)
 - Select “Scan for Devices”
 - It should find your LiveAmp device as long as it is in pairing mode
 - Confirm the connection by searching for “Device Serial Number” on top of the window
 - This SN should match the number on the back of the LiveAmp device
 - Press “Link”
 - You can minimize this window now
- Open Lab Recorder (Shortcut on Desktop)
 - Set your participant number correctly here.
 - We will come back to this window when we are ready to start recording
- Open BrainVision LSL Viewer (Shortcut on Desktop)
 - While this application is not completely necessary in our data collection, it allows the researcher to view the LSL inputs.
- Open Aurora 2023.9 (Shortcut on Desktop)
 - Select “Devices” which is the fifth option across the top (Should be the default open page)
 - Select **both** devices that are connected
 - Select “Select Device” in the bottom right
 - Select “Configurations” which is the last option across the top
 - Select the “Navy” Configuration

- Select “Edit” in the bottom left of the right window showing the Navy configuration.
 - Across the top, the last tab says “LSL Stream Names”
 - Confirm that the info under “Trigger in stream name” says “LSL4Unity.Samples.SimpleCollisionEvent”
 - Press “Save” in the top right
 - Select “Select” in the bottom right of the “Navy” configuration window
- Set Up Quest 3 Headset
 - Confirm that the Oculus App is open
 - Connect the Link Cable from the Oculus to the PC
 - In the headset
 - A window should appear asking if you want to enable or not the quest link
 - Press “Enable”
 - If the window does not appear


After Setting Up Participant with fNIRS / EEG / Cap / Head-Mounted Display (Quest 3)

- Aurora 2023.9
 - Press play to begin testing
- Open BrainVision Recorder (Shortcut - Green on far left)
 - EEG Testing has to occur in this software
- Unity Application
 - Press Play in the top middle



- Aurora 2023.9
 - Across the bottom of the screen, near the bottom middle, a section says “Trigger stream” and the bubble next to it should be green.
 - This indicates that the trigger received from Unity is received
 - If it is still red, then confirm that the play button has been pressed in Unity
 - You can begin recording by pressing the red circle on the top, which is the fifth option now.
- Lab Recorder
 - Confirm Unity and Aurora are both running
 - Press “Update” if you do not see the LSL inputs/outputs as such below
 - You should now see the window below
 - Select the same items shown here
 - “LiveAmpSN-#...”
 - “LiveAmpSN-#...”

- “LSL4Unity.Samples.SimpleCollisionEvent...”
 - Press “Start” to start recording EEG data with LSL inputs being received.



- LSL Viewer
 - Press Connect - should automatically sync and work alongside Lab Recorder
 - Press “Start” if you would like to witness the LSL streams
- Unity Application
 - Continue with the provided instructions (Need to merge documents)

When Complete

- Aurora 2023.9
 - Press the red square, which is the 6th option across the top, to stop recording.
- Lab Recorder
 - Press “Stop” to stop recording
- LSL Viewer

- Press “Stop”
- Unity Application
 - Can press “Start” again to stop the application

Data Collection

- Aurora 2023.9
 - The sixth tab across the top says “Description”
 - At the bottom of this window, you can press the blue link to open the file location in File Explorer.
 - Lsl.tri is the file of the trigger events
 - You can now close the application
- Lab Recorder
 - Note the “Saving to...” location
 - Navigate to this location via the File Explorer
 - Once you have this open in File Explorer, you can close the Lab Recorder
- LSL Viewer
 - No data to be collected from this application.
- Unity Application
 - Unity logs

Data Analysis

- Lsl.tri from File Explorer can currently only be opened from a text editor
 - Does not yet show the metadata - on the to-do list
- MatLab
 - Transfer Files
 - You should already have a file explorer opened that looks similar to this:
“C:\Users\shine\Documents\CurrentStudy\sub-P001\ses-S004\eeg”
 - Open this next to another file explorer opened at
“C:\Users\shine\Documents\MATLAB” to have easy copy and paste
 - Copy the .xdf files from the CurrentStudy folder to the MATLAB location
 - .xdf file example:
“sub-P001_ses-S004_task-Default_run-001_eeg”
 - Allow the Code to Run
 - These files should already be in the MATLAB folder:
 - “Load_xdf_innerloop.mexw64”
 - “Load_xdf.m”
 - If they are not, refer to this GitHub documentation:
<https://github.com/xdf-modules/xdf-Matlab/tree/0cdf054391fff7f0ea3416ee632ad1bd73d6623b>
 - Run the following code:
 - 1st: addpath('C:\Users\shine\Documents\MATLAB')

- 2nd: streams =
load_xdf('sub-P001_ses-S001_task-Default_run-001_eeg.xdf')
- Any additional LSL EEG logs? Be sure to change the name streams to not overwrite previous data
 - For me it was also easy to press the up button to access the previous command and change the participant or session number as needed.
- Actually Analysis (Section Needs Work)
 - Open the saved stream from the workspace on the far right
 - You should see an array of structs
 - Double click one
 - Double click "Info" to see whether the LSL received was from Unity or LiveAmp
 - Double click "Time Series" to see when sent from Unity, the corresponding action
 - For example: "Start Button Pressed"
 - Double click "Time Stamp" to see the corresponding time stamp

Self Documentation - Modify LSL Signals

- VS Code
 - Code for text input:

```
public void lslStuff() {
    if (outlet != null)
    {
        sample[0] = "Start Button Pressed ID:" +
gameObject.GetInstanceID();
        // Debug.Log(sample[0]);
        outlet.push_sample(sample);
    }
}
```

- Continue VS Code
 - Difficult to find, but this could help for future reference if wanting to send different stimulus outputs
 - "ChangePlaneColours"